



LCR Hallcrest Ltd.
Riverside Buildings, Dock Road,
Connah's Quay, Flintshire,
CH5 4DS U.K.
Tel: +44 (0)1244 817107
Fax: +44 (0)1244 818502
Email: sales@lcrhallcrest.com
Web: www.lcrhallcrest.com

Advanced Thermographic Technologies

BriteLine®

Instructions for Use

- Apply Temperature indicating label to outside of urine specimen cup.
- Reversible analog line moves along the scale indicating temperature reading in °C or °F
- Read temperature after line stops moving.