



LCR Hallcrest Ltd.  
Riverside Buildings, Dock Road,  
Connah's Quay, Flintshire,  
CH5 4DS U.K.  
Tel: +44 (0)1244 817107  
Fax: +44 (0)1244 818502  
Email: sales@lcrhallcrest.com  
Web: www.lcrhallcrest.com

---

Advanced Thermographic Technologies

## Brite Event®

---

### Instructions for Use

---

- Apply Temperature indicating label to outside of urine specimen cup.
- Read temperature after colours stop changing
- Green is correct colour
- If Green is not visible, subtract 1° from highest lit event