

BriteLine®

Product Data Sheet

BriteLine® Urine Temperature Indicators

Standard Range: 31.5°C-38.5°C / 89-101°F

Construction:

Adhesive Tape	Adhesive: Single-sided adhesive paper label
Covering Film	Polyester
Color Change Material	Non-toxic mixture of liquid crystal esters
Accuracy	±0.6°C / ±1°F
Shelf Life	12 months from invoice date when stored at room temperature (i.e. 21°C (70°F) & 50% relative humidity).

Application Data: Apply temperature indicating label to outside of urine specimen cup. Reversible analogue line moves along the scale to show detailed temperature reading in °C or °F.

Function: Display of temperature, as indicated by a bright line, of urine specimen cup. Label is adjusted with a 1°F offset.

Use: Ideal for use in clinics, hospitals, and doctors' offices.